

**Introduction to Aseptic Technique (Experiments 6 and 7)**

- Instructor will examine everyone's streak plate and provide suggestions.
- Streak another petri plate (TSA or NA) from your broth culture
- Inoculate another broth culture from your petri plate.
  - Incubate both the petri plate and the broth culture at RT. Appropriately discard your cultures from last week.

**Microscopy (Handout, Experiment 1)**

- Instructor will talk you through the process.

**Microscopy Exercise (Handout)**

- Perform individually. Turn in the completed exercise by the end of class today.