

## Product News

### Cane Attachment

#### *Cane Plus Inc*

Grip 'n Rise, a simple collar attachment for canes, is designed to help patients stand up securely by providing a center of gravity and relieving pressure on the knees. The collar is made from polycarbonate plastic and bears a weight up to 275 lb. It slips onto the cane and locks on with an Allen key.



Circle No 100 on the Reader Service Card or visit [www.apta.org/adinfo](http://www.apta.org/adinfo)

### Run Belt

#### *Fitness With Joe*

The patented EZ Run Belt is designed to help patients poststroke immediately improve mobility, balance and control. The belt gives firm foundation and alignment, says the manufacturer, and it prevents hyperextension of the knee while the patient is walking, standing, or in rehabilitation.



Circle No 101 on the Reader Service Card or visit [www.apta.org/adinfo](http://www.apta.org/adinfo)

*Product information is supplied by manufacturers. Products described should not be construed as having the endorsement of APTA.*

### Portable Biofeedback

#### *Neurogym Technologies Inc*

The Neurogym Trainer is a portable, multifaceted biofeedback



device designed to help patients relearn motor control and coordination. Using inputs from a variety of sensors to play a video game, targeted muscle activity or movement can be encouraged or discouraged.

Circle No 102 on the Reader Service Card or visit [www.apta.org/adinfo](http://www.apta.org/adinfo)

### Wheelchair

#### *Freedom Designs Inc*

The TriPod, a manual wheelchair featuring mid-wheel drive, allows the user's hand to fall naturally into place while propelling the chair. The chair features large wheels and a tilting frame that allows the user or caregiver to maneuver efficiently over small objects and low curbs.



Circle No 103 on the Reader Service Card or visit [www.apta.org/adinfo](http://www.apta.org/adinfo)

### Tilting Table

#### *Tri WG*

Tri WG's Tilt Table is designed specifically for "expanded capacities" and safely accommodates up to 750 lb, with a 90-degree tilt range, 24-inch-wide foot board, and twin locking, heavy-duty casters for easy mobility.



Circle No 104 on the Reader Service Card or visit [www.apta.org/adinfo](http://www.apta.org/adinfo)

### Light Therapy

#### *Dynatronics Corp*

The Dynatron X3TM delivers 3 independent light therapy treatments simultaneously (1 light probe and 2 pads). The physical therapist can choose between 4 Dynatron Solaris Light Probes (combinations of infrared, laser, and blue light) while delivering 2 unattended Dynatron XPTM Light Pad treatments.



Circle No 105 on the Reader Service Card or visit [www.apta.org/adinfo](http://www.apta.org/adinfo)

### Sit-to-Stand Trainer

#### *Neurogym Technologies Inc*

The Sit-to-Stand Trainer enables patients to initiate the movement of standing even with weakened or seemingly unsuccessful effort. Through continued training and gradual reduction of the counterweight, a patient can relearn the "standing-up" skill.



Circle No 106 on the Reader Service Card or visit [www.apta.org/adinfo](http://www.apta.org/adinfo)